Package 'bPeaks'

February 19, 2015

Type Package

Title bPeaks: an intuitive peak-calling strategy to detect

eukaryotic genomes Version 1.2 Date 2014-02-28 Author Jawad MERHEJ and Gaelle LELANDAIS Maintainer Gaelle LELANDAIS <gaelle.lelandais@univ-paris-diderot.fr> Description bPeaks is a simple approach to identify transcription factor binding sites from ChIP-seq data. Our general philosophy is to provide an easy-to-use tool, well-adapted for small eukaryotic genomes (< 20 Mb). bPeaks uses a combination of 4 cutoffs (T1, T2, T3 and T4) to mimic ``good peak" properties as described by biologists who visually inspect the ChIP-seq data on a genome browser. For yeast genomes, bPeaks calculates the proportion of peaks that fall in promoter sequences. These peaks are good candidates as transcription factor binding sites. License GPL Depends R (>= 2.10) NeedsCompilation no Repository CRAN Date/Publication 2014-02-28 17:55:13 R topics documented: bPeaks-package baseLineCalc bPeaksAnalysis dataPDR1 dataReading</gaelle.lelandais@univ-paris-diderot.fr>
Date 2014-02-28 Author Jawad MERHEJ and Gaelle LELANDAIS Maintainer Gaelle LELANDAIS <gaelle.lelandais@univ-paris-diderot.fr> Description bPeaks is a simple approach to identify transcription factor binding sites from ChIP-seq data. Our general philosophy is to provide an easy-to-use tool, well-adapted for small eukaryotic genomes (< 20 Mb). bPeaks uses a combination of 4 cutoffs (T1, T2, T3 and T4) to mimic ``good peak" properties as described by biologists who visually inspect the ChIP-seq data on a genome browser. For yeast genomes, bPeaks calculates the proportion of peaks that fall in promoter sequences. These peaks are good candidates as transcription factor binding sites. License GPL Depends R (>= 2.10) NeedsCompilation no Repository CRAN Date/Publication 2014-02-28 17:55:13 R topics documented: bPeaks-package baseLineCalc bPeaksAnalysis dataPDR1</gaelle.lelandais@univ-paris-diderot.fr>
Author Jawad MERHEJ and Gaelle LELANDAIS Maintainer Gaelle LELANDAIS <gaelle.lelandais@univ-paris-diderot.fr> Description bPeaks is a simple approach to identify transcription factor binding sites from ChIP-seq data. Our general philosophy is to provide an easy-to-use tool, well-adapted for small eukaryotic genomes (< 20 Mb). bPeaks uses a combination of 4 cutoffs (T1, T2, T3 and T4) to mimic ``good peak" properties as described by biologists who visually inspect the ChIP-seq data on a genome browser. For yeast genomes, bPeaks calculates the proportion of peaks that fall in promoter sequences. These peaks are good candidates as transcription factor binding sites. License GPL Depends R (>= 2.10) NeedsCompilation no Repository CRAN Date/Publication 2014-02-28 17:55:13 R topics documented: bPeaks-package baseLineCalc bPeaksAnalysis dataPDR1</gaelle.lelandais@univ-paris-diderot.fr>
Maintainer Gaelle LELANDAIS <gaelle.lelandais@univ-paris-diderot.fr> Description bPeaks is a simple approach to identify transcription factor binding sites from ChIP-seq data. Our general philosophy is to provide an easy-to-use tool, well-adapted for small eukaryotic genomes (< 20 Mb). bPeaks uses a combination of 4 cutoffs (T1, T2, T3 and T4) to mimic ``good peak" properties as described by biologists who visually inspect the ChIP-seq data on a genome browser. For yeast genomes, bPeaks calculates the proportion of peaks that fall in promoter sequences. These peaks are good candidates as transcription factor binding sites. License GPL Depends R (>= 2.10) NeedsCompilation no Repository CRAN Date/Publication 2014-02-28 17:55:13 R topics documented: bPeaks-package baseLineCalc bPeaksAnalysis dataPDR1</gaelle.lelandais@univ-paris-diderot.fr>
Description bPeaks is a simple approach to identify transcription factor binding sites from ChIP-seq data. Our general philosophy is to provide an easy-to-use tool, well-adapted for small eukaryotic genomes (< 20 Mb). bPeaks uses a combination of 4 cutoffs (T1, T2, T3 and T4) to mimic ``good peak" properties as described by biologists who visually inspect the ChIP-seq data on a genome browser. For yeast genomes, bPeaks calculates the proportion of peaks that fall in promoter sequences. These peaks are good candidates as transcription factor binding sites. License GPL Depends R (>= 2.10) NeedsCompilation no Repository CRAN Date/Publication 2014-02-28 17:55:13 R topics documented: bPeaks-package baseLineCalc bPeaksAnalysis dataPDR1
seq data. Our general philosophy is to provide an easy-to-use tool, well-adapted for small eukaryotic genomes (< 20 Mb). bPeaks uses a combination of 4 cutoffs (T1, T2, T3 and T4) to mimic ``good peak" properties as described by biologists who visually inspect the ChIP-seq data on a genome browser. For yeast genomes, bPeaks calculates the proportion of peaks that fall in promoter sequences. These peaks are good candidates as transcription factor binding sites. License GPL Depends R (>= 2.10) NeedsCompilation no Repository CRAN Date/Publication 2014-02-28 17:55:13 R topics documented: bPeaks-package baseLineCalc bPeaksAnalysis dataPDR1
Depends R (>= 2.10) NeedsCompilation no Repository CRAN Date/Publication 2014-02-28 17:55:13 R topics documented: bPeaks-package baseLineCalc bPeaksAnalysis dataPDR1
NeedsCompilation no Repository CRAN Date/Publication 2014-02-28 17:55:13 R topics documented: bPeaks-package baseLineCalc bPeaksAnalysis dataPDR1
Repository CRAN Date/Publication 2014-02-28 17:55:13 R topics documented: bPeaks-package baseLineCalc bPeaksAnalysis dataPDR1
Date/Publication 2014-02-28 17:55:13 R topics documented: bPeaks-package baseLineCalc bPeaksAnalysis dataPDR1
R topics documented: bPeaks-package
R topics documented: bPeaks-package
bPeaks-package
baseLineCalc
dataSmoothing 1 peakDetection 14 peakDrawing 1 peakLocation 19 yeastCDS 20

2 bPeaks-package

Index 22

bPeaks-package bPeaks: an intuitive peak-calling strategy to detect transcription factor binding sites from ChIP-seq data in small eukaryotic genomes

Description

bPeaks is a simple approach to identify transcription factor binding sites from ChIP-seq data. Our general philosophy is to provide an easy-to-use tool, well-adapted for small eukaryotic genomes (< 20 Mb). bPeaks uses a combination of 4 cutoffs (T1, T2, T3 and T4) to mimic "good peak" properties as described by biologists who visually inspect the ChIP-seq data on a genome browser. For yeast genomes, bPeaks calculates the proportion of peaks that fall in promoter sequences. These peaks are good candidates as transcription factor binding sites.

Details

Package: bPeaks Type: Package Version: 1.2

Date: 2014-02-28 License: GPL

Author(s)

Jawad MERHEJ and Gaelle LELANDAIS Maintainer: Gaelle LELANDAIS <gaelle.lelandais@univ-paris-diderot.fr>

References

More information can be found online: http://bpeaks.gene-networks.net/

See Also

```
http://bpeaks.gene-networks.net/
```

```
# get library
library(bPeaks)

# STEP 1: get PDR1 data (ChIP-seq experiments, IP and control samples,
# related to the transcription factor Pdr1 in yeast Saccharomyces
# cerevisiae)
data(dataPDR1)
```

bPeaks-package 3

```
# STEP 2 : bPeaks analysis (only 10 kb of chrIV are analyzed here,
          as an illustration)
bPeaksAnalysis(IPdata = dataPDR1$IPdata[40000:50000,],
               controlData = dataPDR1$controlData[40000:50000,],
      cdsPositions = dataPDR1$cdsPositions,
              windowSize = 150, windowOverlap = 50,
              IPcoeff = 4, controlCoeff = 2,
      log2FC = 1, averageQuantiles = 0.5,
               resultName = "bPeaks_example")
# --> Result files (PDF and BED) are written in the working directory.
# -> bPeaks analysis, all chromosome IV and default parameters (optimized for yeasts)
# STEP 1: get PDR1 data (ChIP-seq experiments, IP and control samples,
# related to the transcription factor Pdr1 in yeast Saccharomyces
# cerevisiae)
data(dataPDR1)
# STEP 2: bPeaks analysis
bPeaksAnalysis(IPdata = dataPDR1$IPdata,
               controlData = dataPDR1$controlData,
      cdsPositions = dataPDR1$cdsPositions,
              windowSize = 150, windowOverlap = 50,
               IPcoeff = 2, controlCoeff = 2,
               log2FC = 2, averageQuantiles = 0.9,
               resultName = "bPeaks_PDR1",
               peakDrawing = TRUE)
# STEP 3 : procedure to locate peaks according to
     gene positions
peakLocation(bedFile = "bPeaks_PDR1_bPeaks_allGenome.bed",
           cdsPositions = yeastCDS$Saccharomyces.cerevisiae,
           outputName = "bPeakLocation_finalPDR1", promSize = 800)
# -> Note that cds (genes) positions are stored in bPeaks package for several yeast
# species
data(yeastCDS)
summary(yeastCDS)
                         Length Class
                                           Mode
# Debaryomyces.hansenii 31370 -none-
                                            character
# Eremothecium.gossypii 23615 -none-
                                            character
# Kluyveromyces.lactis
                         25380 -none-
                                            character
# Pichia.sorbitophila
                          55875 -none-
                                            character
# Saccharomyces.kluyveri 27790 -none-
                                            character
# Yarrowia.lipolytica
                          32235 -none-
                                            character
# Zygosaccharomyces.rouxii 24955 -none-
                                            character
# Saccharomyces.cerevisiae 5 data.frame list
                            5 data.frame list
# Candida.albicans
                             5 data.frame list
# Candida.glabrata
```

4 baseLineCalc

End(Not run)

baseLineCalc

Function to calculate the average number of reads mapped on each nucleotide in a genome

Description

This function calculates the mean genome-wide read depth.

Usage

baseLineCalc(covData)

Arguments

covData

A vector with the numbers of sequences aligned at each genomic position to be considered in the analysis

Details

This function adds the numbers of sequences observed at each position and divides this number by the genome size (total number of nucleotides).

Value

The average number of reads mapped on each nucleotide in the genome.

Note

Detailed information and tutorials can be found online http://bpeaks.gene-networks.net/.

Author(s)

Gaelle LELANDAIS

References

http://bpeaks.gene-networks.net/

See Also

bPeaksAnalysis

bPeaksAnalysis 5

Examples

```
# get library
library(bPeaks)

# get PDR1 data
data(dataPDR1)

# mean genome-wide read depth in IP data
meanIPcov = baseLineCalc(dataPDR1$IPdata[,3])
print(meanIPcov)

# mean genome-wide read depth in control data
meanContCov = baseLineCalc(dataPDR1$controlData[,3])
print(meanContCov)
```

bPeaksAnalysis

Function to run the entire bPeaks procedure

Description

This function allows to detect basic peaks (bPeaks) using the procedure described in the function peakDetection. Chromosomes are analyzed successively. Several values(regarding thresholds T1, T2, T3 and T4 and other parameters) can be specified simultaneously in order to rapidly compare the obtained results and evaluate parameter relevance.

Usage

Arguments

IPdata A dataframe with sequencing results of the IP sample. This dataframe has three

columns (chromosome, position, number of sequences) and should have been

created with the dataReading function

controlData A dataframe with sequencing results of the control sample. This dataframe has

three columns (chromosome, position, number of sequences) and should have

been created with the dataReading function

cdsPositions Not mandatory. A table (matrix) with positions of CDS (genes). Four columns

are required (chromosome, starting position, ending position, strand (W or C), description). CDS positions for several yeast species are stored in bPeaks pack-

age (see the dataset yeastCDS and also peakLocation function)

6 bPeaksAnalysis

smoothingValue The number (n/2) of surrounding positions to use for mean calculation in the

dataSmoothing function

windowSize Size of the sliding windows to scan chromosomes

windowOverlap Size of the overlap between two successive windows

IPcoeff Threshold T1. Value for the multiplicative parameter that will be combined

with the value of the mean genome-wide read depth (see baseLineCalc). As an illustration, if the IPcoeff = 6, it means that to be selected, the IP signal should be GREATER than 6 * (the mean genome-wide read depth). Note that a vector with different values can be specified, the bPeaks analysis will be therefore repeated

using successively each value for peak detection

controlCoeff Threshold T2. Value for the multiplicative parameter that will be combined

with the value of the mean genome-wide read depth (see baseLineCalc). As an illustration, if the controlCoeff = 2, it means that to be selected, the control signal should be LOWER than 2 * (the mean genome-wide read depth). Note that a vector with different values can be specified, the bPeaks analysis will be

therefore repeated using successively each value for peak detection

log2FC Threshold T3. Threshold to consider log2(IP/control) values as sufficiently im-

portant to be interesting. Note that a vector with different values can be specified, the bPeaks analysis will be therefore repeated using successively each value

for peak detection

averageQuantiles

Threshold T4. Threshold to consider $(\log 2(IP) + \log 2(\text{control}))/2$ as sufficiently important to be interesting. This parameter ensures that the analyzed genomic region has enough sequencing coverage to be reliable. These threshold should be between [0, 1] and refers to the quantile value of the global distribution observed

with the analyzed chromosome

resultName Name for output files created during bPeaks procedure

peakDrawing TRUE or FLASE. If TRUE, the function peakDrawing is called and PDF files

with graphical representations of detected peaks are created.

promSize Size of the genomic regions to be considered as "upstream" to the annotated

genomic features (see documentation of the function peakLocation for more

information).

withoutOverlap If TRUE, this option allows to filter peak that are located in a promoter AND a

CDS.

Details

More information together with tutorials can be found online http://bpeaks.gene-networks.net/.

Value

BED files for each chromosomes and a final BED file combining all the results with information regarding detected peaks (genomic positions, mean IP signal, etc.). These files are all saved in the R working directory. Summaries of parameter calculations and peak detection criteria are shown in PDF files (saved in the working directory).

bPeaksAnalysis 7

Note

Detailed information and tutorials can be found online http://bpeaks.gene-networks.net/

Author(s)

Gaelle LELANDAIS

References

```
http://bpeaks.gene-networks.net/
```

See Also

peakDetection dataReading dataSmoothing baseLineCalc peakDrawing peakLocation

```
# get library
library(bPeaks)
# STEP 1: get PDR1 data
data(dataPDR1)
# STEP 2 : bPeaks analysis (only 10 kb of chrIV are analyzed here,
           as an illustration)
bPeaksAnalysis(IPdata = dataPDR1$IPdata[40000:50000,],
               controlData = dataPDR1$controlData[40000:50000,],
               windowSize = 150, windowOverlap = 50,
               IPcoeff = 4, controlCoeff = 2, log2FC = 1,
               averageQuantiles = 0.5,
               resultName = "bPeaks_example",
               peakDrawing = TRUE, promSize = 800)
## Not run:
# STEP 2 : bPeaks analysis (all chromosome)
bPeaksAnalysis(IPdata = dataPDR1$IPdata, controlData = dataPDR1$controlData,
                cdsPositions = dataPDR1$cdsPositions,
                smoothingValue = c(20),
                windowSize = c(150), windowOverlap = 50,
                IPcoeff = c(2), controlCoeff = c(2), log2FC = c(2),
                averageQuantiles = c(0.9),
                resultName = "bPeaks_PDR1_chr4",
                peakDrawing = TRUE, promSize = 800)
# To repeat the bPeaks analysis with different parameters
bPeaksAnalysis(IPdata = dataPDR1$IPdata, controlData = dataPDR1$controlData,
                cdsPositions = dataPDR1$cdsPositions,
                smoothingValue = c(20),
                windowSize = c(150), windowOverlap = 50,
                IPcoeff = c(2, 4, 6), controlCoeff = c(2, 4, 6), log2FC = c(2, 3),
                averageQuantiles = c(0.7, 0.9),
                resultName = "bPeaks_PDR1_chr4_paremeterEval",
```

8 dataPDR1

```
peakDrawing = FALSE, promSize = 800)
 # -> Summary table is created and saved as "peakStats.Robject" in the working directory
 # as well as a text file named "_bPeaks_parameterSummary.txt"...
 load("peakStats.Robject")
 # This table comprises different information regarding peak detection (number of peaks,
 # mean size of peaks, mean IP signal, mean control signal, etc.)
 peakStats[1:2,]
        smoothingValue windowSize windowOverlap IPcoeff controlCoeff log2FC
 # [1,]
                     20
                               150
                                               50
                                                        1
                                                                             1
 # [2,]
                     20
                               150
                                               50
                                                                      1
                                                        1
        averageQuantiles bPeakNumber meanSize meanIPsignal meanControlSignal
 # [1,]
                                  308 209.091
                                                     276.047
                                                                         71.534
                      0.5
                                                     287.808
                                                                         74.002
 # [2,]
                      0.7
                                  294 205.782
       meanLog2FC bPeakNumber_beforeFeatures bPeakNumber_afterFeatures
 #
 # [1,]
              1.571
                                             99
                                                                        80
 # [2,]
              1.589
                                             94
                                                                        77
       bPeakNumber_inFeatures
 # [1,]
 # [2,]
                             53
 ## End(Not run)
dataPDR1
                          ChIP-seq results (IP and control samples) obtained with the transcrip-
```

Description

ChIP-seq experiments were performed in order to identify the genomic regions that interact with the transcription factor Pdr1, in yeast Saccharomyces cerevisiae. Two samples (IP and control) were sequenced simultaneously using the Illumina technology (ENS transcriptome platform). Only the data for chrIV are available here, but complete datasets can be found online: http://bpeaks.gene-networks.net

tion factor Pdr1 in yeast Saccharomyces cerevisiae

Usage

data(dataPDR1)

Format

dataPDR1\$IPdata: IPdata and controlData are dataframes with three columns. The first column comprises chromosome names, the second column comprises base positions and the third column comprises the numbers of sequences mapped at the considered position. dataPDR1\$controlData: IPdata and controlData are dataframes with three columns. The first column comprises chromosome names, the second column comprises base positions and the third column comprises the numbers of sequences mapped at the considered position. dataPDR1\$cdsPositions: A table with annotated positions of genes in yeast S. cerevisiae. The first column indicates chromosome names, the

dataPDR1

second and third columns indicate respectively "start" and "end" positions of genes, and the fourth column indicates the gene annotation (according to the Saccharomyces Genome Database (SGD http://www.yeastgenome.org/)).

Details

Complete procedure to analyze sequencing data (intial FASTQ files) can be found online:

http://bpeaks.gene-networks.net. Initial read length was 50 bases. After quality controls and filtering of low quality bases, around 30.000.000 of sequences (IP sample) and around 88.000.000 of sequences (control sample) were independently mapped on the genome using the bowtie algorithm [1]. Output files (SAM format) were converted into BAM files and indexed using the SAMTOOLS suite [2]. Numbers of sequences mapped on each nucleotide in the reference genome were finally calculated using the "genomeCoverageBed" tool available from the BEDTOOLS suite [3].

Source

Sample sequencing was performed at the "transcriptome platform", ENS institute in Paris (France), http://www.transcriptome.ens.fr.

References

More information concerning this dataset can found online: http://bpeaks.gene-networks.net.

- [1] Langmead B, Trapnell C, Pop M, Salzberg SL. Ultrafast and memory-efficient alignment of short DNA sequences to the human genome. Genome Biol 10:R25.
- [2] Li H.*, Handsaker B.*, Wysoker A., Fennell T., Ruan J., Homer N., Marth G., Abecasis G., Durbin R. and 1000 Genome Project Data Processing Subgroup (2009) The Sequence alignment/map (SAM) format and SAMtools. Bioinformatics, 25, 2078 9. [PMID: 19505943]
- [3] Quinlan AR and Hall IM, 2010. BEDTools: a flexible suite of utilities for comparing genomic features. Bioinformatics. 26, 6, pp. 841 842.

```
# get library
library(bPeaks)
# get data
data(dataPDR1)
summary(dataPDR1)
             Length Class
# IPdata
              3 data.frame list
# controlData 3
                     data.frame list
# cdsPositions 5
                     data.frame list
# run peak calling, comparing IP and control samples
# (only 10 kb of chrIV are analyzed here, as an illustration)
bPeaksAnalysis(IPdata = dataPDR1$IPdata[40000:50000,],
              controlData = dataPDR1$controlData[40000:50000,],
              windowSize = 150, windowOverlap = 50,
              IPcoeff = 4, controlCoeff = 2,
```

10 dataPDR1

```
log2FC = 1, averageQuantiles = 0.5,
              resultName = "bPeaks_example",
              peakDrawing = TRUE, promSize = 800)
## Not run:
# -> bPeaks analysis, all chromosome IV and default parameters (optimized for yeasts)
# STEP 1: get PDR1 data (ChIP-seq experiments, IP and control samples,
# related to the transcription factor Pdr1 in yeast Saccharomyces
# cerevisiae)
data(dataPDR1)
# STEP 2: bPeaks analysis
bPeaksAnalysis(IPdata = dataPDR1$IPdata,
              controlData = dataPDR1$controlData,
      cdsPositions = dataPDR1$cdsPositions,
              windowSize = 150, windowOverlap = 50,
              IPcoeff = 2, controlCoeff = 2,
              log2FC = 2, averageQuantiles = 0.9,
              resultName = "bPeaks_PDR1",
              peakDrawing = TRUE)
# STEP 3 : procedure to locate peaks according to
     gene positions
peakLocation(bedFile = "bPeaks_PDR1_bPeaks_allGenome.bed",
           cdsPositions = yeastCDS$Saccharomyces.cerevisiae,
           outputName = "bPeakLocation_finalPDR1", promSize = 800)
# -> Note that cds (genes) positions are stored in bPeaks package for several yeast
# species
data(yeastCDS)
summary(yeastCDS)
                         Length Class
                                           Mode
# Debaryomyces.hansenii 31370 -none-
                                            character
# Eremothecium.gossypii 23615 -none-
                                            character
# Kluyveromyces.lactis
                         25380 -none-
                                            character
                          55875 -none-
# Pichia.sorbitophila
                                            character
# Saccharomyces.kluyveri 27790 -none-
                                            character
# Yarrowia.lipolytica
                          32235 -none-
                                            character
# Zygosaccharomyces.rouxii 24955 -none-
                                          character
# Saccharomyces.cerevisiae 5 data.frame list
# Candida.albicans
                             5 data.frame list
                            5 data.frame list
# Candida.glabrata
# number of detected peaks
print(resLocation$numPeaks)
# number of peaks "upstream" annotated genes (W strand)
print(resLocation$beforeFeatures)
# number of peaks "in" annotated genes
print(resLocation$inFeatures)
```

11 dataReading

```
# number of peaks "upstream" annotated genes (C strand)
print(resLocation$afterFeatures)
## End(Not run)
```

dataReading

Function to import sequencing results in R

Description

Sequencing results should be converted to datafiles with numbers of sequences mapped on each nucleotide in the reference genome. These files can be generated from indexed BAM files (mapping results) using the "genomeCoverageBed" tool available from the BEDTOOLS suite [1]. More information concerning file descriptions can be found online:

http://bpeaks.gene-networks.net/.

File format should be as follows (chromosome, position, number of sequences):

Usage

dataReading(IPfile, controlFile, yeastSpecies = NULL)

Arguments

IPfile Name of the file with sequencing results related to IP sample controlFile Name of the file with sequencing results related to control sample yeastSpecies

Not mandatory. Annotations to be used for locations of peaks in promoters. An-

notations of CDS are available in bPeaks for yeasts: Debaryomyces.hansenii,

Eremothecium.gossypii, Kluyveromyces.lactis, Pichia.sorbitophila, Saccharomyces.kluyveri,

Yarrowia.lipolytica, Zygosaccharomyces.rouxii, Saccharomyces.cerevisiae, Can-

dida.albicans, Candida.glabrata (see data yeastCDS)

Details

To obtain a required file from a BAM file (resultFile.bam), the command line is (SHELL):

12 dataReading

```
genomeCoverageBed -ibam resultFile.bam -d > resultFile.txt

More information concerning file conversions can be found online:

http://bpeaks.gene-networks.net/.
```

Value

A list with three elements (\$IPdata, \$controlData, \$cdsPositions): IP data, control data and (if specified by user) CDS positions for locations of peaks in promoters.

Note

Conversion of file formats regarding sequencing results can be a tricky task. Detailed information can be found online http://bpeaks.gene-networks.net/. Don't hesitate to contact us for further discussions.

Author(s)

Gaelle LELANDAIS

References

[1] Quinlan AR and Hall IM, 2010. BEDTools: a flexible suite of utilities for comparing genomic features. Bioinformatics. 26, 6, pp. 841 842.

See Also

peakLocation dataPDR1 yeastCDS

```
# get library
library(bPeaks)
## Not run:
# Sequencing result files associated to PDR1 datasets (IP and control files)
# can be downloaded from our website http://bpeaks.gene-networks.net/.
# They are respectively named "IP_genomeCoverage.txt" (IP sample) and
# "INPUT_genomeCoverage.txt" (control sample).
# Import in R the sequencing results using S. cerevisiae CDS annotations.
data(yeastCDS)
seqResult = dataReading("IP_genomeCoverage.txt", "INPUT_genomeCoverage.txt",
yeastSpecies = yeastCDS$Saccharomyces.cerevisiae)
# IP data
seqResult$IPdata
# control data
seqResult$controlData
# run peak detection from IP and control data (with default parameters)
```

dataSmoothing 13

dataSmoothing

Function to smooth sequencing coverage along a chromosome

Description

This function allows to obtain a smoothed signal of the genome-wide read depth. Simple moving average (SMA) procedure is used. At each genomic position, the sequencing coverage is replaced by the unweighted mean of the n surrounding positions (n/2 before and n/2 after).

Usage

```
dataSmoothing(vecData, widthValue = 20)
```

Arguments

vecData A vector with the numbers of sequences aligned at each genomic position to be

considered in the analysis

width Value The number (n/2) of surrounding positions to use for mean calculation

Details

Detailed information and tutorials can be found online http://bpeaks.gene-networks.net/.

Value

A vector with the smoothed signal. Note that the SMA procedure creates missing values at the beginning and at the end of the vector with the smoothed signal.

Note

Detailed information and tutorials can be found online http://bpeaks.gene-networks.net/.

Author(s)

Gaelle LELANDAIS

References

```
http://bpeaks.gene-networks.net/
```

See Also

bPeaksAnalysis baseLineCalc

14 peakDetection

Examples

```
# get data
data(dataPDR1)
# inital IP signal
iniIPsignal = dataPDR1$IPdata[,3]
par(mfrow = c(2,2))
# plot initial IP signal
plot(iniIPsignal[416900:417400], type = "h",
    xlab = "genomic position", ylab = "sequencing coverage",
   main = "IP sample (PDR1 data)\nno smoothing",
   col = "red")
# calculate and plot smoothed signal (widthValue = 5)
smoothedIPsignal = dataSmoothing(vecData = iniIPsignal, widthValue = 5)
plot(smoothedIPsignal[416900:417400], type = "h",
    xlab = "genomic position", ylab = "sequencing coverage",
   main = "IP sample (PDR1 data)\nsmoothing (widthValue = 5)",
   col = "pink")
# calculate and plot smoothed signal (widthValue = 10)
smoothedIPsignal = dataSmoothing(vecData = iniIPsignal, widthValue = 10)
plot(smoothedIPsignal[416900:417400], type = "h",
    xlab = "genomic position", ylab = "sequencing coverage",
   main = "IP sample (PDR1 data)\nsmoothing (widthValue = 10)",
   col = "pink")
# calculate and plot smoothed signal (widthValue = 20)
smoothedIPsignal = dataSmoothing(vecData = iniIPsignal, widthValue = 20)
plot(smoothedIPsignal[416900:417400], type = "h",
   xlab = "genomic position", ylab = "sequencing coverage",
   main = "IP sample (PDR1 data)\nsmoothing (widthValue = 20)",
   col = "pink")
```

peakDetection

Peak calling method, i. e. identification of genomic regions with a high density of sequences (reads)

Description

bPeaks uses a sliding window to scan the genomic sequence. Four criterion define interesting regions: 1) a high number of reads in the IP sample (T1 = IPthreshold); 2) a low number of reads in the control sample (T2 = controlThreshold); 3) a high value of log(IP/control) (T3 = ratioThreshold) and 4) a good sequencing coverage (IP and control samples) (T4 = averageThreshold). Parameters for peak detection, therefore, rely on four threshold values (T1, T2, T3 and T4). Graphical outputs and BED files are provided allowing the user to rapidly assess the relevance of the chosen parameters.

peakDetection 15

Usage

Arguments

IPdata A dataframe with sequencing results of the IP sample. This dataframe has three

columns (chromosome, position, number of sequences) and should have been

created with the dataReading function

controlData A dataframe with sequencing results of the control sample. This dataframe has

three columns (chromosome, position, number of sequences) and should have

been created with the dataReading function

chrName Name of the chromosome to be scanned with the sliding windows (to compare

IP and control signals and detect interesting regions)

windowSize Size of the sliding window to scan chromosomes windowOverlap Size of the overlap between two successive windows

outputName Name for output files created during bPeaks procedure

baseLineIP Value of the mean genome-wide read depth (IP sample). This value is calculated

using the baseLineCalc function

baseLineControl

Value of the mean genome-wide read depth (control sample). This value is

calculated using the baseLineCalc function

IPthreshold Threshold T1. Threshold to consider IP signal as sufficiently important to be

interesting. Note that these threshold is a multiplicative parameter that will be combined with the calculated "baseLineIP" (see before) value. As an illustration, if the IPthreshold = 6, it means that to be selected, the IP signal should be

GREATER than 6 * baseLineIP

controlThreshold

Threshold T2. Threshold to consider control signal as sufficiently low to be interesting. Note that these threshold is a multiplicative parameter that will be combined with the calculated "baseLineControl" (see before) value. As an illustration, if the controlThreshold = 2, it means that to be selected, the control

signal should remain LOWER than 2 * baseLineControl

ratioThreshold Threshold T3. Threshold to consider log2(IP/control) values as sufficiently important to be interesting

averageThreshold

Threshold T4. Threshold to consider $(\log 2(IP) + \log 2(\text{control}))/2$ as sufficiently important to be interesting. These parameter is important to ensure that the analyzed genomic region has enough sequencing coverage to be reliable. These threshold should be between [0, 1] and refers to the quantile value of the global distribution observed with the analyzed chromosome

peakDrawing TRUE or FLASE. If TRUE, the function peakDrawing is called and PDF files with graphical representations of detected peaks are created.

16 peakDetection

Details

Detailed description of the bPeaks procedure together with tutorials can be found online: http://bpeaks.gene-networks.net/.

Value

A matrix with genomic positions of the detected peaks. Summaries of parameter calculations and peak detection criteria are shown in PDF files (saved in the working directory).

Note

Detailed information and tutorials can be found online http://bpeaks.gene-networks.net/. Don't hesitate to contact us for further discussions.

Author(s)

Gaelle LELANDAIS

References

```
http://bpeaks.gene-networks.net/
```

See Also

dataReading baseLineCalc peakDrawing bPeaksAnalysis

```
# get library
library(bPeaks)
# get PDR1 data
data(dataPDR1)
# combine IP and control data
allData = cbind(dataPDR1$IPdata, dataPDR1$controlData)
colnames(allData) = c("chr", "pos", "IPsignal", "chr", "pos", "controlSignal")
# calculate baseline IP and control values
        = baseLineCalc(allData$IPsignal)
print(paste("Baseline coverage value in IP sample : ", round(lineIP, 3)))
lineControl = baseLineCalc(allData$controlSignal)
print(paste("Baseline coverage value in control sample : ", round(lineControl, 3)))
print("")
# get list of chromosomes
chromNames = unique(allData[,1])
# start peak detection on the first chromosome
```

peakDrawing 17

```
print(paste("Starting analysis of chromosome ", chromNames[1]))
# information for one chromosome
subData = allData[allData[,1] == chromNames[1],]
# only 10 kb are analyzed here (as an illustration)
         = subData[40000:50000,3]
vecControl = subData[40000:50000,6]
# smooth of the data
smoothedIP
           = dataSmoothing(vecData = vecIP, widthValue = 20)
smoothedControl = dataSmoothing(vecData = vecControl, widthValue = 20)
# peak detection
detectedPeaks = peakDetection(IPdata = smoothedIP, controlData = smoothedControl,
               chrName = as.character(chromNames[1]),
               windowSize = 150, windowOverlap = 50,
               outputName = paste("bPeaks_example_", chromNames[1], sep = ""),
               baseLineIP = lineIP, baseLineControl = lineControl,
               IPthreshold = 4, controlThreshold = 2,
               ratioThreshold = 1, averageThreshold = 0.5,
               peakDrawing = TRUE)
# print detected genomic positions
print(detectedPeaks)
```

peakDrawing

Function to draw graphical representations of genomic regions detected using bPeaks methodology

Description

This function allows to create PDF files, with graphical representations of the detected basic peaks (bPeaks). Genomic regions are shown together with the values of the parameters used to detect the region.

Usage

Arguments

vecIP	Vector with sequencing depth at each nucleotide (from start pos = 1 to end)
vecControl	Vector with sequencing depth at each nucleotide (from start pos = 1 to end)
lineIP	Threshold values used for peak detection (IP signal)
lineControl	Threshold values used for peak detection (control signal)

18 peakDrawing

lineFC	Threshold val	ues used for peak	detection (log2)	(IP/control) values)

lineAverage Threshold values used for peak detection (average log2(IP) and log2(control)

values)

posInf Genomic position to start the representation posSup Genomic position to end the representation

add Number of bases before and after posInf and posSup to add

title Graphic main title

Details

More information can be found online: http://bpeaks.gene-networks.net/.

Value

Image in x11() terminal

Note

Detailed information and tutorials can be found online http://bpeaks.gene-networks.net/.

Author(s)

Gaelle LELANDAIS

References

```
http://bpeaks.gene-networks.net/
```

See Also

bPeaksAnalysis

peakLocation 19

peakLocation	Function to locate detected basic peaks (bPeaks) according to predifined chromosomal features

Description

Starting from a BED file with positions of detected peaks and a table with positions of CDS (genes), this function allows to dentify the peaks that are located "upstream" or "in" annotated CDS. Annotations of CDS for different yeast species are available in bPeaks package (see data yeastCDS).

Usage

Arguments

bedFile	Name of a BED file with positions of detected peaks (using bPeaks or another program)
cdsPositions	A table (matrix) with positions of CDS (genes). Four columns are required (chromosome, starting position, ending position, strand (W or C), description)
withoutOverlap	If TRUE, this option allows to filter peak that are located in a promoter AND a CDS.
outputName	Name for output files
promSize	Genomic size to be considered as promoter (upstream to CDS)

Details

More information can be found online http://bpeaks.gene-networks.net/.

Value

Graphics and text files (saved in the R working directory).

Note

Detailed information and tutorials can be found online http://bpeaks.gene-networks.net/.

Author(s)

Gaelle LELANDAIS

References

```
http://bpeaks.gene-networks.net/
```

20 yeastCDS

See Also

bPeaksAnalysis dataReading yeastCDS

Examples

```
## Not run:
# -> bPeaks analysis with (all chromosome and default parameters optimized for yeasts)
# STEP 1: get PDR1 data and annotations in yeasts
data(dataPDR1)
data(yeastCDS)
# STEP 2: bPeaks analysis
bPeaksAnalysis(IPdata = dataPDR1$IPdata,
               controlData = dataPDR1$controlData,
               windowSize = 150, windowOverlap = 50,
               IPcoeff = 6, controlCoeff = 4,
               log2FC = 2, averageQuantiles = 0.9,
               resultName = "bPeaks_PDR1",
               peakDrawing = TRUE)
# STEP 3 : procedure to locate peaks according to
# predefined chromosomal features
peakLocation(bedFile = "bPeaks_PDR1_bPeaks_allGenome.bed",
            cdsPositions = yeastCDS$Saccharomyces.cerevisiae,
   withoutOverlap = FALSE,
            outputName = "bPeakLocation_finalPDR1", promSize = 800)
## End(Not run)
```

yeastCDS

Annotations of CDS for different yeast species

Description

Annotations of CDS for different yeast species. Data were collected (July 2013) from: Genolevures http://genolevures.org/, SGD http://www.yeastgenome.org/, CGD http://www.candidagenome.org/databases.

Usage

```
data(yeastCDS)
```

Format

yeastCDS\$Saccharomyces.cerevisiae: A table (matrix) with positions of CDS (genes) in yeast S. cerevisiae. Four columns are required (chromosome, starting position, ending position, strand (W or C), description).

yeastCDS 21

Details

This R object is a list composed of multiple tables (matrices) - one for each species - with positions of CDS (genes).

Source

Data were collected (July 2013) from Genolevures http://genolevures.org/, SGD http://www.yeastgenome.org/and CGD http://www.candidagenome.org/databases.

References

```
Genolevures database http://genolevures.org/
SGD database http://www.yeastgenome.org/
CGD database http://www.candidagenome.org/
```

```
# get library
library(bPeaks)
# get data
data(yeastCDS)
# different species for wich information is available
summary(yeastCDS)
                                        Mode
                        Length Class
# Debaryomyces.hansenii 31370 -none-
                                       character
# Eremothecium.gossypii 23615 -none-
                                          character
# Kluyveromyces.lactis 25380 -none-
                                          character
# Pichia.sorbitophila 55875 -none-
                                          character
# Saccharomyces.kluyveri 27790 -none-
                                          character
# Yarrowia.lipolytica 32235 -none-
                                          character
# Zygosaccharomyces.rouxii 24955 -none-
                                          character
# Saccharomyces.cerevisiae 5 data.frame list
                            5 data.frame list
# Candida.albicans
                           5 data.frame list
# Candida.glabrata
# CDS annotations for yeast Debaryomyces hansenii
yeastCDS$Debaryomyces.hansenii[1:10,]
      chrm start end
                              strand geneName
# [1,] "Deha2A" "2023" "6370" "C"
                                     "DEHA2A00110g"
# [2,] "Deha2A" "6587" "7810" "C"
                                     "DEHA2A00132g"
# [3,] "Deha2A" "8314" "9354" "W"
                                     "DEHA2A00154g"
# [4,] "Deha2A" "9632" "9844" "C"
                                     "DEHA2A00176g"
# [5,] "Deha2A" "13806" "14132" "W"
                                     "DEHA2A00198g"
# [6,] "Deha2A" "14558" "16519" "C"
                                     "DEHA2A00220g"
# [7,] "Deha2A" "17520" "19442" "W"
                                     "DEHA2A00242g"
# [8,] "Deha2A" "22619" "23977" "C"
                                     "DEHA2A00264g"
# [9,] "Deha2A" "24949" "25434" "W"
                                     "DEHA2A00286g"
#[10,] "Deha2A" "26440" "26640" "C"
                                     "DEHA2A00308g"
```

Index

*Topic ChIP-seq results dataPDR1, 8 yeastCDS, 20	*Topic signal smoothing dataSmoothing, 13 *Topic small eukaryotic genome
*Topic ChIP-seq	bPeaks-package, 2
bPeaks-package, 2 bPeaksAnalysis, 5 peakDetection, 14 *Topic Pdr1 transcription factor	baseLineCalc, 4, 6, 7, 13, 15, 16 bPeaks (bPeaks-package), 2 bPeaks-package, 2 bPeaksAnalysis, 4, 5, 13, 16, 18, 20
dataPDR1, 8 yeastCDS, 20 *Topic Saccharomyces cerevisiae dataPDR1, 8 yeastCDS, 20	dataPDR1, 8, <i>12</i> dataReading, <i>5</i> , <i>7</i> , 11, <i>15</i> , <i>16</i> , <i>20</i> dataSmoothing, <i>6</i> , <i>7</i> , 13
*Topic bPeaks peakDrawing, 17 *Topic data reading dataReading, 11	peakDetection, 5, 7, 14 peakDrawing, 6, 7, 15, 16, 17 peakLocation, 5-7, 12, 19
*Topic deep sequencing bPeaks-package, 2 bPeaksAnalysis, 5 peakDetection, 14	yeastCDS, 5, 11, 12, 19, 20, 20
*Topic genome coverage baseLineCalc, 4 *Topic moving average	
dataSmoothing, 13 *Topic peak calling	
<pre>bPeaks-package, 2 bPeaksAnalysis, 5 peakDetection, 14</pre>	
*Topic protein binding sites bPeaks-package, 2 bPeaksAnalysis, 5	
peakDetection, 14 *Topic protein-DNA interactions	
bPeaks-package, 2 bPeaksAnalysis, 5 peakDetection, 14 *Topic read depth	
hacelineCale 4	